



Technical Bulletin

Ref: **TB-14006**

Issue date: **February 2022**

Product: INOmax DS_{IR}[®]

Priority: Low

Affected parts: Injector Module

Classification: Information for
Distribution

Subject: **Part Number Change**

WARNING: Use only parts/accessories designated for use with the INOmax DS_{IR}.

This technical bulletin is written to advise clinicians that the injector module part number for the INOmax DS_{IR} has changed.

The injector module has been modified to comply with RoHS requirements (hazardous substance free). The current injector module will still function with the INOmax DS_{IR}, but providing new injector modules outside of the United States will require the RoHS compliant injector module beginning July 22, 2014.

The outward appearance will not change other than the part number on the packaging.

Part/Accessory	Part Number
Injector Module, RoHS Compliant	90726

If you have further questions, please contact technical support at 1-877-566-9466 (North America) or your specific country manager.

Mallinckrodt Manufacturing LLC
6603 Femrite Drive,
Madison, WI 53718-6801 USA

Mallinckrodt, the "M" brand mark and the Mallinckrodt Pharmaceuticals logo are trademark of a Mallinckrodt company. Other brands are trademarks of a Mallinckrodt company or their respective owner. ©2016 Mallinckrodt



Europe

European Authorized Representative
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

